


Technical Data Sheet

Rev 01 09/2024

TF.01.11

1



CRS® Custom Ortho Model Resin

Printer compatibility 385 nm / 405 nm

CHARACTERISTICS

PROPERTY	VALUE	UNIT
Colour	Light Grey	-
Viscosity	200-250	mPa.s @ 25 °C
Density	1.11-1.14	g/cm3

PERFORMANCE FEATURES

Property	Result	Unit	Standard
Flexural Strength	78	MPa	ASTM D790
Flexural Modulus	2365	MPa	ASTM D790
Tensile Strength	35.3	MPa	ASTM D638
Tensile Modulus	2005	MPa	ASTM D638
Shore Hardness	83	Shore D	ASTM D2240

No color change is observed after curing.